




INFUSE YOUR PETS WITH A WORK OF HEART



PETite

Vet Infusion Pump with Fluid Warmer



Exclusive to




Pumps



With convenience in mind, select from 3 versatile positions:






Flexible Setting Adjustments:

Temperatures can be continuously adjusted between 86°F and 104°F, with a control accuracy of $\pm 1.8^\circ\text{F}$.




Multiple Safety Features:

For safe and secure infusion process, the heat automatically shuts off when the temperature exceeds 122°F.



Rapid Heating:

Equipped with four powerful heat units, this pump rapidly achieves the desired temperature.



Extended Heating Tube

The 1100mm-long heating tube efficiently warms the medicinal liquid while minimizing the loss of drug activity.


Safe and Reliable

Designed with a flexible three-chamber structure, this unit is fully sealed and features a waterproof design.



Portable Installation:

Includes a heater mounting base for convenient storage and organization.



Precise Temperature Control:

Equipped with dual NTC temperature sensors for precise and reliable heating control.

