

ABC Model 4100 Syringe Infusion Pump

Always Best CureTM

Ease of Operation Particularly for **Enteral** Feedings

User Customization promotes friendly programming

Multiple Enteral Delivery Modes

- Various Delivery modes (mL/hr, Volume/Time)
- Adjustable Flow Rate ranges
- Accuracy +/- 3% excluding syringe variations

Several Enteral Safety Messages

- ENTERAL Message is displayed on LCD when an enteral manufacturer is selected
- User must respond to specific ENTERAL question(s)

Accommodate Various Enteral Syringes

Manufacturers with Automatic Size Sensing

- Enteral – 1 ml to 100 ml

Easy Syringe Loading

- Front Syringe Loading with Syringe Light
- Easy Plunger Loading with Plunger Sensor

Multiple Sensors to Monitor Delivery Accuracy

- Syringe Size
- Syringe Plunger Force
- Syringe Plunger Movement
- Motor Shaft Rotation

Metal Housing Design

- Non-plastic housing to reduce breakage
- Improved EMI immunity (IEC 60601-1-2)
- Splash Proof for more effective and easier cleaning

Large Graphic LCD Display with White Backlight (promotes easy viewing)

Syringe Light can be viewed from afar (delineates infusion status)

Unique Labeling area can be used to indicate infusion information (see example above)

Orange Keypad color is optional (see above)

AC / Battery operated device with Externally Removable Battery Pack

- Standard Battery Capacity (12 hours at 5 ml/hr)
- No External Power Adaptor required

LED & LCD Alerts with Adjustable Audio Alarm Volume and Frequency

RS232 Interface Capabilities

Integrated Carrying Handle that can attach to a 360° Rotating Universal Pole or Rail Mount

041612

